e-sportshop.cz

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
MANUÁL

ID: 65409

EAN: 8591792654097

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Quadriceps Femoris		Circle muscle	30s-45s
Rectus femoris	20s-45sx2 times	Straight knee position	15s-30s
Lateral femur	20s-45sx2 times	Flexed knee position	15s-30s
Medial femur	20s-45sx2 times	Calf triceps	
Middle femoral muscle	20s-45sx2 times	Overhead muscles	20s-30sx3s
		Halibut	20s-30sx3times

The direction of movement
is from top to bottom and
the patient slowly bends
the knee slightly

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Small round head vertical down

Lateral humeral ankle tendon 30S40S

If you are afraid of hitting the bone surface, you can squeeze the muscle with
your hand and hit it downward

Small round head downward vertically:

Flexor carpi radialis	10s-15sx3 times
Palmaris longus	12s-15sx3 times
Flexor carpi radialis ulnaris	10s-15sx3 times
Finger extensors	5s-10sx2 times
Flexor carpi radialis short	5s-10sx2 times
Flexor carpi radialis longus	5s-10sx2 times

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Triceps	15s-20sx3reps
Elbow muscles	10s-20sx2 sets
Superior trapezius bundle	15s-20sx3 reps
Middle fascicle	15s-20sx3 reps
Inferior fasciculus	15s-20sx3 reps
Latissimus dorsi	30s-1minx2 reps
Direction from the middle to the sides, from top to bottom	
Small round head vertical down	
Supraspinatus	10s-15sx3reps
Scapular lift	10s-12sx2 times
Large circular muscle	10s-12sx2 reps
Small round muscle	10s-12sx2 times
Infraspinatus	10s-12sx2 times
Direction from the middle to the sides	

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scle (middle and lower fibers)

Middle fasciculus	15s-20sx3 reps
Inferior fasciculus	15s-20sx3 reps
Rhomboids	20s-30sx3 reps
Vertical spinal muscles (lumbar segment)	30s-1minx4 reps
latissimus dorsi	30s-1minx4
Thoracolumbar fascia	1minx3 times
Sacral joint	1minx3 times
Small flat head vertical down	
Truncalis major	30s-1minx6 times
Supine position	
Internal and external obliques	10s-15sx2 times
Small round head vertical down	
erector spinae (thoracic segment)	15s-20sx2 times
Lumbar square	20s-3.0sx2 times
Iliopsoas muscle	30s-40sx4 times

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s maximus	30s-1minx4 times
Gluteus minimus and skeleton	30s-1minx3 times
Sacral skeletal joint	30s-1minx2 times
Lateral position, bend the straightened leg below by 90	
Small flat head down vertically:	
Gluteus medius and skeleton	30s-1r
Sacral skeletal joint	30s-1r
Skeletal tibial bundle	30s-1r
Around greater trochanter	30s-1r
Supine position (unilateral)	
Skeletal muscles	30s-1minx3 times
Slow flexion of the patient's hip	
Skeletal margin-anterior skeleton	30s-1minx2 times
Rectus femoris inguinalis	1minx3 times

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Direction: inward upward to outward downward, along the muscle

Rousseau muscle	
Biceps femoris	20s-30sx3 times
Semitendinosus	20s-30sx3 times
Semimembranosus	20s-30sx3 times
From top to bottom	

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Orbital cord muscle	
Biceps femoris	20s-30sx3 reps
Semitendinosus	220s-30sx3 times
Semimembranosus muscle	220s-30sx3 times
Direction of movement from top to bottom	
U-shaped head vertical downward circling muscle	
Straight knee position	215s-30s
Flexed knee position	215s-30s
Calf triceps	
Gastrocnemius	220s-30sx3 reps
Hallux valgus	220s-30sx3 reps
Four-point position with forefoot on the ground	
Small round head perpendicular to the sole of the foot	
Plantar fascia	220s-30sx3 times

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